

## Using an Osstell® together with a Multipeg®

Measured ISQ	Add	Measured ISQ	Add	Measured ISQ	Add
20	-4	44	0	68	-1
21	-4	45	-1	69	0
22	-4	46	-1	70	0
23	-4	47	-1	71	0
24	-4	48	-1	72	1
25	-3	49	-2	73	1
26	-3	50	-2	74	2
27	-3	51	-2	75	2
28	-3	52	-2	76	2
29	-3	53	-2	77	3
30	-3	54	-2	78	3
31	-2	55	-3	79	4
32	-2	56	-3	80	4
33	-2	57	-3	81	4
34	-2	58	-3	82	5
35	-2	59	-2	83	5
36	-2	60	-2	84	5
37	-1	61	-2	85	6
38	-1	62	-2	86	6
39	-1	63	-2	87	6
40	-1	64	-2	88	7
41	-1	65	-2	89	7
42	-1	66	-1	90	8
43	-1	67	-1		

Take the value in the "Add" column and add to the measured ISQ

